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" Driving up to Mendocino and thinking about an experiment to look at one particular letter of the genetic code, I designed a system in my mind. As I repaired the things I thought could go wrong with it, suddenly I generated something that if I did it over and over again would be PCR. It would go 2, 4, 8, 16, 32 . . . in 30 cycles make as many base pairs from one little region as I had in the whole genome! That was the eureka point. I said holy shit! By putting the triphosphates [DNA building blocks] in there myself, I could do this process over and over and amplify the DNA.

I slammed on the brakes and stopped by the side of the road to calculate it out. Jennifer objected groggily to the delay and the light, but I exclaimed I had discovered something fantastic. Unimpressed, she went back to sleep.



A couple of miles down the road I stopped again. I realized I could use these bastards, the oligonucleotides [short pieces of DNA], and get the enzymes to reproduce as big a piece as I wanted to. They didn't have to be aimed at just one base pair. Hell, I could do a whole sequence. I realized you can cut the sequence out from a great big molecule. Pretty cool! Just cut, paste, and amplify.»

Cetus Corporation sold the patent for PCR to Hoffman-LaRoche for the staggering \$300 million - the most money ever paid for a patent. Mullis meanwhile received a \$10,000 bonus.

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



Thermus aquaticus

Yellow stone National Park, California

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<http://biotools.nubic.northwestern.edu/OligoCalc.html>

Primer-BLAST

A tool for finding specific primers

Finding primers specific to your PCR template (using Primer3 and BLAST).

Primers for target on one template

[Primers common for a group of sequences](#)

[Reset page](#)

[Save search parameters](#)

[Retrieve recent results](#)

[Publication](#)

[Tips for finding specific primers](#)

PCR Template

Enter accession, gi, or FASTA sequence (A refseq record is preferred) [Clear](#)

```
>gi|34328764:16431-17671 Mytilus galloprovincialis haplotype M mitochondrion, rRNA
till tRNA
AAAGCATTAACTATAAAAGATACCCCTAAATGTTTTGTAAAAAGCGGGCTAAACAAGGATAGTATAAT
AAATTTTAGTACCTTTTGATAAGGGTTTTCAAGACAAATTTAAGTATTTTCCCGAATGAAAGA
GAGTTATTTGTAGAGTTAAAAATCGTGGTAAAGATTTAATTAATTATAAAATAGCGGCTAGATACTAT
```

Or, upload FASTA file

Choose File

No file chosen

Range [Clear](#)

From

To

Forward primer

Reverse primer

Primer Parameters

Use my own forward primer (5'->3' on plus strand)

[Clear](#)

Use my own reverse primer (5'->3' on minus strand)

[Clear](#)

PCR product size

Min

Max

of primers to return

Primer melting temperatures (T_m)

Min

Opt

Max

Max T_m difference

Exon/intron selection

A refseq mRNA sequence as PCR template input is required for options in the section [Clear](#)

Exon junction span

Exon junction match

Min 5' match

Min 3' match

Max 3' match

Note: Parameter values that differ from the default are highlighted in yellow

Primer Pair Specificity Checking Parameters

Specificity check

Enable search for primer pairs specific to the intended PCR template [?](#)

Search mode

Automatic [?](#)

Database

Custom [?](#)

Enter accession number, gi, or FASTA sequence [Clear](#)

AY363687.2

Or, upload file: No file chosen

Exclusion

Exclude predicted Refseq transcripts (accession with XM, XR prefix) Exclude uncultured/environmental sample sequences [?](#)

Organism

Mytilus galloprovincialis (taxid:29158)

Enter an organism name (or organism group name such as enterobacteriaceae, rodents), taxonomy id or select from the suggestion list as you type. [?](#)

[Add more organisms](#)

Entrez query (optional)

[?](#)

Primer specificity stringency

Primer must have at least total mismatches to unintended targets, including

at least mismatches within the last bps at the 3' end. [?](#)

Ignore targets that have or more mismatches to the primer. [?](#)

Max target size

[?](#)

Allow splice variants

Allow primer to amplify mRNA splice variants (requires refseq mRNA sequence as PCR template input) [?](#)

Show results in a new window Use new graphic view [?](#)

Oligo Calc: Oligonucleotide Properties Calculator

Enter Oligonucleotide Sequence Below
OD calculations are for single-stranded DNA or RNA

Nucleotide base codes

TGC CCT AGA GGC GTA GAA GCC TC

Reverse Complement Strand(5' to 3') is:

GAG GCT TCT ACG CCT CTA GGG CA

5' modification (if any)

3' modification (if any)

Select molecule

nM Primer

Measured Absorbance at 260 nanometers

mM Salt (Na⁺)

Calculate

Swap Strands

BLAST

mfold

Physical Constants

Length: Molecular Weight: GC content: %

1 ml of a sol'n with an Absorbance of at 260 nm

is microMolar and contains micrograms.

Melting Temperature (T_M) Calculations

°C (Basic)

°C (Salt Adjusted)

°C (Nearest Neighbor)

Thermodynamic Constants Conditions: 1 M NaCl at 25°C at pH 7.

RlnK cal/(°K*mol)

deltaH Kcal/mol

deltaG Kcal/mol

deltaS cal/(°K*mol)

Deprecated Hairpin/self dimerization calculations

(Minimum base pairs required for single primer self-dimerization)

(Minimum base pairs required for a hairpin)

Check Self-Complementarity

CR-F: 5'-TGCCCTAGAGGCGTAGAAGCCTC -3' (23 nt)

Minimum base pairs required for single primer self-dimerization: 5.
Minimum base pairs required for a hairpin: 4.

Potential hairpin formation :

None !

3' Complementarity:
None !

All potential self-annealing sites are marked in red (allowing 1 mis-match):

None !

Oligo Calc: Oligonucleotide Properties Calculator

Enter Oligonucleotide Sequence Below
OD calculations are for single-stranded DNA or RNA

Nucleotide base codes

CTC TGA GAA ATG CTT ATC AGC TG

Reverse Complement Strand(5' to 3') is:

CAG CTG ATA AGC ATT TCT CAG AG

5' modification (if any)

3' modification (if any)

Select molecule

ssDNA

50 nM Primer

50 mM Salt (Na⁺)

1 Measured Absorbance at 260 nanometers

Calculate

Swap Strands

BLAST

mfold

Physical Constants

Length: Molecular Weight: GC content: %

1 ml of a sol'n with an Absorbance of at 260 nm

is microMolar and contains micrograms.

Melting Temperature (T_M) Calculations

1 °C (Basic)

2 °C (Salt Adjusted)

3 °C (Nearest Neighbor)

Thermodynamic Constants Conditions: 1 M NaCl at 25°C at pH 7.

RlnK cal/(°K**mol*)

deltaH Kcal/*mol*

deltaG Kcal/*mol*

deltaS cal/(°K**mol*)

Deprecated Hairpin/self dimerization calculations

(Minimum base pairs required for single primer self-dimerization)

(Minimum base pairs required for a hairpin)

Check Self-Complementarity

Citation: Kibbe WA. 'OligoCalc: an online oligonucleotide properties calculator'. (2007)
Nucleic Acids Res. **35**(webserver issue): May 25. ([Abstract/Full text](#))

This page may be freely linked or distributed for any educational or non-commercial use.

CR-R: 5'-CTCTGAGAAATGCTTATCAGCTG-3' (23 nt)

Minimum base pairs required for single primer self-dimerization: 5.
Minimum base pairs required for a hairpin: 4.

Potential hairpin formation :

5' CTCTGAGAAATGCTTATCAGCTG 3'

3' Complementarity:
None !

All potential self-annealing sites are marked in red (allowing 1 mis-match):

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